Kat no: BS-NA-513-100

vNAT® Transfer Tüpü

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Package Insert

1. Product Content and Shelf Life

Product	Content	Definition	Amount
Bio-Speedy® vNAT® Transfer Tube*	Bio-Speedy® vNAT® Transfer Tube	Microbial nucleic acid storage and stabilization device	100 tubes

*Shelf life: 12 months

2. Storage and Transport Conditions of Bio-Speedy® vNAT® Transfer Tube

Bio-Speedy® vNAT® Transfer Tube should be stored and transported at **15-50°C**.

3. Intended Use

Bio-Speedy® vNAT® Transfer Tube, 2 mL of viral nucleic acid extractive and preservative liquid. When clinical specimens suspected of respiratory tract infection are transferred in Bio-Speedy® vNAT® Transfer Tube, the liquid inside the tube can be used directly in real-time PCR (qPCR) reactions. The nucleic acid extractive and preservative liquid inactivates all viral, bacterial, or eukaryotic pathogens in the sample, 5 minutes after contact with the clinical specimen.

WARNING: On the web page linked with the QR code, Instruction for Use (IFU) is given. Do not start the analysis without reading the IFU.

4. Sampling Protocol

Swab samples:

Place the *Specimen Collection Swab* into the Bio-Speedy® vNAT® Transfer Tube after sampling.

Samples in liquid form:

1. After vortexing the sample tube for 15 seconds, transfer liquid sample (1mL ≤ sample volume ≤ 2 mL) into the Bio-Speedy® vNAT® Transfer Tube.

Sputum:

2. Equal amounts of sputum and nuclease-free water (VTM or 1xPBS can be used instead of water) are mixed and homogenized. It is kept for 1 minute to settle the solid/ viscous layer. Transfer the sample from the upper liquid part (1mL ≤ sample volume ≤ 2 mL) into the Bio-Speedy® vNAT® Transfer Tube.

5. Transportation, Storage and Analysis Protocol

1. Bio-Speedy® vNAT® Transfer Tube containing the sample should be delivered and stored to the laboratory at **2-8°C**.
2. Vortex the Bio-Speedy® vNAT® Transfer Tube at the highest speed for 5 seconds.
3. The sample in the Bio-Speedy® vNAT® Transfer Tube is ready to use in qPCR.

6. Manufacturer and Technical Support

Bioeksen R&D Technologies Incorporated Company

Address: ITU ARI Teknokent 3, 4/B-105, 34467. Sariyer, Istanbul, TURKEY.

Phone: +90 (212) 285 10 17, **Fax:** +90 (212) 285 10 18

Web: www.bioeksen.com.tr, **e-mail:** info@bioeksen.com.tr, support@bioeksen.com.tr